

# GeLC-MS/MS

Standardised nomenclature for gel slice labels in GeLC experiments:

**MAKE SURE THAT LABELLING ON TUBES AND GEL SCHEME ARE THE SAME!**

|          | SAMPLE | SAMPLE | SAMPLE | SAMPLE | SAMPLE | SAMPLE | SAMPLE | SAMPLE | SAMPLE |        |
|----------|--------|--------|--------|--------|--------|--------|--------|--------|--------|--------|
|          | A      | B      | C      | D      | E      | F      | G      | H      | I      |        |
| (MW)     | A1     | B1     | C1     | D1     | E1     | F1     | G1     | H1     | I1     | top    |
| (MARKER) | A2     | B2     | C2     | D2     | E2     | F2     | G2     | H2     | I2     |        |
|          | A3     | B3     | C3     | D3     | E3     | F3     | G3     | H3     | I3     |        |
|          | A4     | B4     | C4     | D4     | E4     | F4     | G4     | H4     | I4     |        |
|          | A5     | B5     | C5     | D5     | E5     | F5     | G5     | H5     | I5     | bottom |

